



***Evanthia:***  
***Leading Filler Technology***

Aesthetic Medicine made in Italy

Evanthia from Bioformula combines advanced technology with the safest manufacturing process - made in Italy, resulting in a safe and effective cross-linked hyaluronic acid dermal filler line for the face, that provides instant, long lasting and natural looking results.

## ***Evanthia 15***

Softens the appearance of fine lines.

Suggested treatments:

- Small periocular wrinkles
- Small perioral wrinkles



## ***Evanthia 20***

Redefines the face without losing its natural expressions.

Suggested treatments:

- Lip contour & enhancement
- Nasolabial folds
- Eyebrow lift
- Facial reshaping
- Marionette lines



## ***Evanthia 25***

Enhances the face volume naturally and non-invasively.

Suggested treatments:

- Volumetric facial reshaping
- Lip enhancement
- Prejowl sulcus
- Deep nasolabial folds
- Temples
- Cheek augmentation
- Chin augmentation
- Rhinofiller
- Jawline contouring





We do not compromise when it comes to quality. Every stage of the production process is carefully monitored to ensure all parameters of the final product meet the conformity limits, with some being way below the accepted limits.\*

- Cross-linking agent residual:  
< 0,05 ppm
- Endotoxin content:  
 $\leq 0,5$  Eu/ml
- Residual protein concentration:  
0,1%

\* Conformity parameters  
set by FDA and EU Regulation 2017/745.



Certifications of our manufacturing partner:  
UNI EN ISO 9001:2015 since 1995 with certificate No. 423  
UNI CEI EN ISO 13485:2016 since 2008 with certificate No. 13790  
ISO 22716:2007 since 2013 with certificate No. 18549



*Evanthia*  
*Suggested*  
*Treatments*



# Evanthia

## Technical Specifications\*

\* Please refer to the I.F.U. for full specifications.

\*\* The estimated average duration depends on various factors: patient's skin type, depth of wrinkles, injection area and the injected quantity. It is recommended to have the first two sessions within 4-6 months to achieve optimal results.

	<i>Evanthia 15 Basic</i>	<i>Evanthia 20 Medium</i>	<i>Evanthia 25 Superior</i>
Product code	40166	40167	40168
Composition	Reticulated hyaluronic acid of non-animal origin	Reticulated hyaluronic acid of non-animal origin	Reticulated hyaluronic acid of non-animal origin
Cross-linking agent	BDDE	BDDE	BDDE
HA Concentration	<u>15 mg/ml</u>	<u>20 mg/ml</u>	<u>25 mg/ml</u>
Raw material molecular weight	2400 - 4000 kDa	2400 - 4000 kDa	2400 - 4000 kDa
Packaging	2 x 1 ml	2 x 1 ml	2 x 1 ml
Blister	2 pre-filled, graduated, disposable, sterile glass syringes with Luer-Lok™ adapter  4 sterile disposable needles 30G 1/2" (CE 0123)	2 pre-filled, graduated, disposable, sterile glass syringes with Luer-Lok™ adapter  4 sterile disposable needles 27G 1/2" (CE 0123)	2 pre-filled, graduated, disposable, sterile glass syringes with Luer-Lok™ adapter  4 sterile disposable needles 27G 1/2" (CE 0123)
Injection depth	Superficial derma	Medium derma	Deep derma
Estimated duration of results**	Up to 6 months	Up to 8 months	Up to 10 months



Bioformula injectable solutions are  
Class III Medical Devices (CE 0373).  
For medical use only.

Bioformula S.r.l. Viale delle Industrie, 8  
20865 Usmate-Velate (MB) Italy +39 039 631 751

